



Floating Ball Blanket for Water Baths

Polycarbonate balls are an ideal way to contain heat and minimise evaporation for a bath when it is impractical to apply a lid.

Model	Description
FBB-200	Bag of 200 Balls suits 10 litre size bath
FBB-400	Bag of 400 balls suits 20 litre size bath



WiseBath Water Bath

(WiseBath) Digital Fuzzy Control water bath with LCD Display and Jog Shuttle Switch for quick temperature adjustments. Stainless steel bath for superior durability and high thermal efficiency. Digital timer function (99hr 59min) with delay time and operation time. Alarm function to indicate errors and timer end.

Models	WB-6	WB-11	WB-22
Bath Dims	300x155x150	302x240x150	500x295x150
Overall Dims	338x195x260	338x280x260	536x335x260
Capacity	6 litre	11 litre	22 litre
Temp Range	Ambient + 5°C to 100°C		
Controller	Digital Fuzzy Controller		
Safety	Adjustable Over Temp / Sensor Error Detect / Low Level Detect		
Options	Lift off S/S gable lid		



Stainless Steel Gable Lid

\* Optional Stainless Steel Gable Lid available for water baths.

- WB-6-GL : Lid for 6 litre bath
- WB-11-GL : Lid for 11 litre bath
- WB-22-GL : Lid for 22 litre bath

**Baths, Circulated**



WiseCircu Circulating Baths

(WiseCircu) Digital Fuzzy Control circulating bath with LCD Display and Jog Shuttle Switch for quick temperature adjustments. Powerful circulation pump to ensure high temperature uniformity. Pump out facility to allow for connection with external instruments. Temperature compensation function supported. Safety devices include over temperature, over current and low water warning.

Models	WCB-6	WCB-11	WCB-22
Bath Dims (wxdxh)	300x155x150	302x240x150	500x295x150
Overall Dims (wxdxh)	338x195x400	338x280x400	536x335x400
Capacity	6 litre	11 litre	22 litre
Temp Range	Ambient + 5°C to 120°C		
Controller	Digital Fuzzy Control, System		
Safety	Adjustable Over Temp / Sensor Error Detect / Low Level Detect		
Options	Lift off S/S gable lid		
Heater	800W	1200W	2000W
Pump	14L/Min (internal/external)		
Timer	99hr 59min		

\* Optional Stainless Steel Gable Lid available for water baths.

- WCB-6-GL : Lid for 6 litre bath
- WCB-11-GL : Lid for 11 litre bath
- WCB-22-GL : Lid for 22 litre bath