



ALUF-370

This unit offers 370 litres of freezer space operating down to  $-22^{\circ}\text{C}$  with manual defrost eliminating fluctuating temperature. Ideal for laboratory applications where a constant temperature is essential. This unit has 6 fixed wire shelves. Standard model has digital temperature control and audible hi/lo temperature alarms.

External finish is in a white PVC coated steel and the interior is plastic. Door has a key lock.

Model	ALUF-370
Capacity	370 litres
Temp. Range	$-18^{\circ}\text{C}$ to $-22^{\circ}\text{C}$
Storage Containers	-
Ext. Dims. (hwxwd mm)	1840x595x600
Int. Dims. (hwxwd mm)	1530x475x440
External Material	White enamel
Internal Material	White Plastic
Electrical	240V 50Hz
Weight	65kg
Options	External connection for BMS



TLUF-520-2D

This double door freezer can operate down to  $-22^{\circ}\text{C}$  with manual defrost eliminating fluctuating temperature. Ideal for laboratory applications where a constant temperature is essential. The two door design ensures that the low temperatures are contained while retrieving items from either the top or bottom of the freezer.

Standard model has digital temperature control and audible hi/lo temperature alarms. External and internal finish is in stainless steel.

Door lock can be fitted to the top door only allowing for the safe storage of valuable items. Additional shelving also available.

Model	TLUF-520-2D
Capacity	520 litres
Temp. Range	$-18^{\circ}\text{C}$ to $-22^{\circ}\text{C}$
Shelves	3 adjustable
Ext. Dims. (hwxwd mm)	1830x640x800
Int. Dims. (hwxwd mm)	2 sections 690x540x700
External Material	Stainless Steel
Internal Material	Stainless Steel
Electrical	240V 50Hz
Weight	120kg
Options	External connection for BMS, additional shelving, top door lock