



205-RM



210-RM

Simultaneous rotating and rocking action provides fast and gentle mixing. All rollers are made of special plastic for non-slip and for easy cleaning. The quiet motor, durable plastic gears, smooth mixing action result in a low noise level. Baseplate is made of stainless steel for easy maintenance.

Model	205RM-5	210RM-10
Roller	5 each	10 each
Speed	40 rpm	50 rpm
Roller Length	325 mm	325mm
Dimension (lxwxh mm)	480x220x110	480x420x110
Net Weight	4.5 kg	8.5 kg



KJMR-II

Kang Jian blood mixer provides simultaneous rotating and rocking for fast and gentle mixing. All 7 rollers are made of plastic for non-slip and for easy cleaning. The quiet motor, durable plastic gears, smooth mixing action result in a low noise level.

Model	KJMR-II
Roller	7 each
Speed	40 rpm (Adjustable)
Roller Length	240 mm
Dimension (lxwxh mm)	420x220x110
Net Weight	4 kg

Mixers, Vortex



Vortex Mixers
250VM / 260VM

Both models have variable speed control which allows slow speed shaking action up to high speed vortexing. Both models have two modes of operation; 'touch mode' by pressing the head of the unit and 'full mode' for continuous operation.

Model 250VM has a single cup head and is suitable for mixing the content of the test tubes.

Model 260VM has 90mm dia flat disc style head suitable for mixing the content of beakers or flasks.

Net weight: 3kg
Speed: 0~3000 rpm
Power: 240V/50Hz