



ALR-850

These five models of Refrigerators are suitable for general laboratory products needed to be stored at a temperature between 2°C and 8°C. All models have forced air design and are fitted with digital controls with hi/low temperature alarms. Low temp safety thermostats can be fitted to prevent freezing in the case of malfunction. Door locks can also be added. Shelves are PVC coated, open wire type and are adjustable. Exterior finish is powder coated steel and interior is white PVC.

Model	ALR-210	ALR-358	ALR-506	ALR-660	ALR-850
Ext Dims (WxDxH cm)	51x56x176	61x64x194	61x74x207	76x74x207	126x72x208
Int Dims (WxDxH cm)	40x43x114	50x53x130	50x59x143	66x60x143	114x59x125
Temp. Range	2°C to 8°C				
Temp. Precision	±0.1°C				
Door	Glass	Glass	Glass	Glass	Glass
Volume (litres)	210	358	506	660	850
Shelf	5	5	5	5	10
Weight (kg)	78	100	118	142	180