



AV061P

The Kaltis range of ultra low temperature cabinets are specifically designed to meet the essential reliability demand of ultra low temperature freezers while ensuring long term durability.

- These freezers are simple to use and maintain and are very quiet while operating.
- CFC free insulation.
- Heated door seal and lid gaskets form multi-point seals, minimizing frost build up especially around access areas.
- Upright and chest cabinets use a down feed evaporator coil which completely surrounds the inner chamber. The coil is positioned for uninterrupted refrigerant flow and optimal cooling efficiency.
- The air cooled refrigeration system has been designed to maintain the set temperature even in elevated ambient temperature environments.
- Key pad, multifunctional display, microprocessor control system.
- Audible alarms alert users in the event of power failure, temperature changes and door ajar.
- Heavy gauge steel cabinet with high impact powder coated exterior finish.
- Stainless steel interior.
- Rounded interior finish for easy cleaning.
- Complete inventory system including stainless steel racks (horizontal for uprights and vertical for chests) and cardboard plus plastic boxes available.
- Optional Kaltis circular chart recorder
- Optional CO2 backup system
- 2 year compressor warranty

Upright Freezer Models

Model	AV039P	AV050P	AV061P	AV073P
Capacity	390 Litres	499 litres	616 litres	732 litres
Temp. Range	-80°C			
Interior Dims	60x50x129	60x50x129	74x64.5x129	88x64.5x129
Exterior Dims	85x74x200	85x88.5x200	99x88.5x200	113x88.5x200
Weight	300	330	350	370
Electrical	240V 50Hz			

Freezer Accessories

Part Number	Description
3121C-01	Kaltis 7 day temperature recorder (150mm dia) chart factory fitted
KRCP	Spare Charts, Box of 50