



HF-151-UV



HF-212-UV

This range of CO₂ Incubators are ideal for cell and tissue culture work and feature:

- Attractive powder coated exterior
- Auto start function
- Microprocessor control system
- Independent overtemp protection
- Accurate temperature control $\pm 0.1^{\circ}\text{C}$ from ambient plus 5°C to $+50^{\circ}\text{C}$ range.
- Accurate CO₂ sensor 0 to 20% range
- Air jacket heating for all over even heating of controlled space.
- High humidity levels maintained
- Stainless steel liner has radiussed corners for easy cleaning
- UV sterilisation (30 minute cycle)
- Rapid recovery time of conditions
- Inner glass door and heated outer door

ORDER CODE	DIMS INT (HxWxD)	LITRES	DIMS EXT (HxWxD)	SHELVES	ELECTRICAL (220V 50Hz)	WEIGHT (kg)
HF-151-UV	60x47x53	151	87x62x68	3	600W	90
HF-212-UV	60x60x59	212	80x91x76	3	770W	125
ADDITIONAL EXTRAS						
HF-151-S	Additional shelf for model HF-151-UV					
HF-151-UVL	Spare UV lamp for model HF-151-UV					
HF-212-S	Additional shelf for model HF-212-UV					
HF-212-UVL	Spare UV lamp for model HF-212-UV					
HF-EA	Electronic CO ₂ Analyser					