



WOV-30

Vacuum Drying Ovens, common applications include suitability for Drying, Baking, Conditioning, Curing, Out-Gassing of solids and liquids, Vacuum embedding, moisture testing and age testing. Available in 30 litre and 70 litre volumes. All controls conveniently situated in the front panel for easy access. Robust and reliable design these vacuum ovens will offer many years of service.

- Electronic digital control technology for accurate control
- User select auto or manual run applications
- Large Vacuum gauge for easy viewing
- High temperature rating silicon door gasket
- Safety glass view window
- Storage ability for alarms, timer on digital control panel
- Lock mode
- Stainless steel 304 internal sealed shell
- Removable inner shelves

Safety features include over temperature protection, over current protection and sensor fail indication.

Model	WOV-30	WOV-70
Capacity	30 litres	70 litres
Internal Size (wxdxh mm)	310x310x310	410x410x410
External Size (wxdxh mm)	680x450x500	780x560x590
Heater Capacity	600W x 2	700W x 2
Type	Natural Convection	
Temp. Range	Ambient to 200°C	
Temp. Accuracy	±2°C	
Temp. Controller	Digital Fuzzy Controller	
Vacuum Range	0 to 75cmHg	
Timer Range	99Hr to 59Minutes	
Timer Controller	Digital LCD Display	
Material Internal	Stainless Steel	
Material External	Powder coated steel	
Insulation	Ceramic Wool	
Door Packing	High Temperature Silicone Rubber	
Shelves	Chrome plated steel	
Window	High Temperature Safety Glass	
Safety Device	Over temperature, over current protection	
Power	AC 230V, 50Hz	