



RK-1D

High precision digital speed control
 Programmable operation modes supported 10 steps , 6 programs
 Alarm for error or timer end
 Convenient setting of timer and speed using jog shuttle switch
 Digital LCD display for timer, RPM and power bar graph
 Storage function for RPM and Timer (99hr 59 min)
 Locking mode
 Smooth start and change of speed
 Over load protection

Model	RK-1D	RK-2D
Platform Dims (wxd mm)	290x200	300x300
Overall Dims (wxdxh mm)	290x240x150	300x330x150
Tilt Angle	9° (default)	
Motor	35W geared DC motor	
Speed Range	5 - 50 digital 100 Hr digital	
Timer	100 Hr digital	
Display	Digital LCD display	
Safety	Over load detector	
Options	Optional platform with 4 pillar	
Power	220V 50Hz	
Alarms	Emergency alarm (buzzer)	

Shakers, Twisting



VS-96TW

Dimpled rubber top allows different sized tubes to be used simultaneously.
 Gentle motion ensures noise and vibration-free operation.

Model	VS-96TW
Speed Range	6 ~ 64 oscillations per min
Speed Controller	electronic controller
Drive motor	brushless geared motor
Angle	0 ~ 45°
Motion	Orbital Rocking Motion
Tray Size	260x235mm
Dims	260Wx240Dx250H mm mm
Mat	Dimpled Rubber Mat
Power	220V 50Hz