



Analogue Soxhlet Water Baths available with 4 or 6 holes with concentric rings. Each bath is constructed of stainless steel and features support clamps for larger flasks. Digital bath control is optionally available.

| Model                    | SWB-4            | SWB-6          |
|--------------------------|------------------|----------------|
| Type                     | 4 hole opening   | 6 hole opening |
| Capacity                 | 11 litre         | 14 litre       |
| Bath Size (wxdxh mm)     | 330x320x100      | 890x150x100    |
| External Size (wxdxh mm) | 365x340x290      | 920x180x190    |
| Heater Capacity          | 1000W            | 1500W          |
| Temp Range               | Ambient to 100°C |                |

Soxhlet Water Baths