


Cold Plates (-5 °C to -15°C)

AUSTRALIAN MADE



FEATURES

- Manufactured by Thermoline our Cold Plates are designed specifically for the Pathology and Histology applications.
- However this by no means limits there usage, if you need to cool, freeze, or maintain that sample below zero for short periods of time you might be surprised how easy our Cold Plates complete your tasks.
- Three models to meet your needs 2 having step down plates to accommodate user being seated at the work bench
- Aluminium work platform
- Large Plate area
- Fast cool down rates
- Thermostat temperature control between -5°C and -15°C
- CFC free refrigeration gas

Safety Devices
√ Over Current Protection

ORDER CODE	DIM. EXT LxWxH cm	PLATE AREA cm	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg.
TCP-1	39x37x20	36x15	1.5A	15
TCP-2	43x35x29	38x30	3A	24
TCP-3	43x68x25	38x30	3A	28