

Freezers, Economy (-15°C to -20°C)

AUSTRALIAN MADE



TEUF-181



TEUF-90

FEATURES

Thermoline economy style upright freezers to operate at -20°C with manual defrost. Ideal for low temperature storage of specimens that cannot tolerate a rise in temperature for defrost. Storage containers of clear plastic. Standard model has digital control and audible hi/lo temperature alarms.

- Manual defrost ensures low temperatures are always maintained during operation
- Fully adjustable open wire shelving (additional shelving available)
- Key lock fitted to door
- Added benefit of digital control with indication of internal conditions
- Easy installation, simply plug in power and set temperature.

Safety Features

- ✓ Alarms fitted with rear contacts
- ✓ Over current protection

Optional Features

TCR-115, 7 day circular chart recorder

TS-1000, upgraded alarms system with built in TEST feature, battery back up and power loss notification

TPD-1, Automatic phone dialer to call when alarms are activated (even your mobile phone)

*POWER GUARD, offering protection against Low Voltage, High Voltage Spikes & Limited Lighting *NEW PRODUCT**

Model	TEUF-90	TEUF-181
Capacity	90 Litres	181 litres
Temp. Range	-15°C to -20°C	
Storage Containers	4 clear plastic	6 clear plastic
Ext. Dims. (hwxwd mm)	870x540x600	1440x545x570
Int. Dims. (hwxwd mm)	650 / 430x400x400	1240x405x380
External Material	White enamel	
Internal Material	White PS	
Electrical	240v 50Hz, Max 4A	
Weight	40kg	60 kg
Options	External connection for BMS (in-house monitoring). Key locks for doors.	