

Freezers (-15°C to -20°C)

AUSTRALIAN MADE



TUF-800-20-2-GD

(shown with optional TS-1000 alarm system and recorder)

* Shown with optional TCR-115 and TS-2000

UPRIGHT STYLE FAN FORCED

SD= Solid Door
GD= Glass Door



TUF-530-20-1-SD
TUF-530-20-1-GD



TUF-800-20-2-SD
TUF-800-20-2-GD

FEATURES

Thermoline have designed two -20°C premium freezers available in solid or double glazed door, suited for long or short term storage of your products and samples. Features include.

- Auto-Defrost feature means no more messy defrosting during critical times
- Select between 530Lt and 800Lt upright freezers (TUF-800-series has 2 doors)
- Fully adjustable open wire shelving (additional shelving available)
- Self closing access door/s
- Fan forced air distribution ensures even temperatures throughout the internal volume
- Added benefit of digital control with indication of internal conditions
- Alarms fitted for added security includes audible and visual signals with mutable function, also supplied with remote contacts at the rear for connections to BMS system (alarm probe thermally lagged to replicate product temperature and reduce annoying false alarms)
- Easy installation, with swivel & lockable castors, simply plug in power and set temperature.

Safety Features

- ✓ Alarms fitted with rear contacts
- ✓ Over current protection
- ✓ Lockable castors

Optional Features

TCR-115, 7 day circular chart recorder
 TS-1000, upgraded alarms system with built in TEST feature, battery back up and power loss notification
 TPD-1, Automatic phone dialer to call when alarms are activated {even your mobile phone}
 POWER GUARD, offering protection against Low Voltage, High Voltage Spikes & Limited Lighting *NEW PRODUCT*

ORDER CODE	DIM. INT HxWxD cm	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg.
TUF-530-20-1-SD	146x58.5x56	530	200x74x71	4	4A	150
TUF-530-20-1-GD	146x58.5x56	530	200x74x71	4	4A	150
TUF-800-20-2-SD	146x85x56	800	200x100x71	4	5A	190
TUF-800-20-2-GD	146x85x56	800	200x100x71	4	5A	190

Upright Freezers (-20°C to -35°C)

AUSTRALIAN MADE



TUF-350-35-1-SD

FEATURES

Thermoline have produced this freezer for products that need to be continuously stored below -25°C, for example Plasma or Blood products. The available size is 350Lt in an upright configuration.

- Hermetically sealed refrigeration system for extremely quiet and reliable operations
- Ample storage capacity of 350Lt
- Wide temperature control range -20°C to -35°C
- Fixed plate shelving, this means your product is in closest contact possible with the freezing coils at all times
- Digital control, featuring digital indication of internal conditions
- Fully Insulated cabinet construction for energy efficiency
- Designed with compressor and electronics in the top compartment above the storage area, this means the heat generated by the system does not impede on internal temperatures.
- Hi/Lo temperature alarms fitted as standard for added product security

Safety Features

- ✓ Alarms fitted with rear contacts
- ✓ Over current protection
- ✓ Lockable castors

Optional Features

TCR-115, 7 day circular chart recorder

TS-1000, upgraded alarms system with built in TEST feature, battery back up and power loss notification

TPD-1, Automatic phone dialer to call when alarms are activated (even your mobile phone)

*POWER GUARD, offering protection against Low Voltage, High Voltage Spikes *NEW PRODUCT**

ORDER CODE	DIM. INT HxWxD cm	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg
TUF-350-35-1-SD	111x55x58	350	173x71x78	3	5A	130