


## Economy Incubators (Ambient +5°C to 60°C)

AUSTRALIAN MADE



### FEATURES

Thermoline's Economy Incubators are designed and manufactured in Australia by our in house engineers. We have 4 models for you to select, making them an excellent economical choice in all non critical laboratory incubation application. Thermoline also manufactures a wide range of fan forced incubators should your needs be for tighter control conditions and larger volumes (see page I4).

#### Standard Features includes

- Stainless steel internal liner (easy clean and non corrosive)
- Removable open wire chrome plated shelves
- Magnetic door gasket means positive door closure
- Light weight and portable
- Digital Microprocessor control complete with dual digital display
- Control error indication
- Simple to use
- Unique "Thermowall" heating system with electronic switching
- Control range ambient +5°C to +60°C
- Control Accuracy at set point  $\pm 0.25^\circ\text{C}$  or better

#### Safety Devices

- ✓ Electrical over current protect
- ✓ Control error codes

#### Optional Accessories

*IGD- Internal glass door installed (prevents loss of internal conditions when access door is opened)*  
*GVW- Glass View Window installed in access door*

ORDER CODE	DIM. INT HxWxD cm	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg	OPTIONAL INNER GLASS DOOR	OPTIONAL GLASS VIEW WINDOW
TEI-13G	22x25x24	13	37x30x30	2	60W	8	IGD	GVW
TEI-43G	50x29x29	43	65x34x35	4	120W	16	IGD	GVW
TEI-70G	60x29x40	70	75x34x46	5	180W	24	IGD	GVW
TEI-100G	65x39x40	100	83x46x47	6	240W	32	IGD	GVW