

Laboratory Refrigerators, Economy (+4°C to +10°C)

AUSTRALIAN MADE




TELR-130-1-SD

TELR-500-1-GD

UPRIGHT STYLE TEMPERATURE RANGE +4°C TO +10°C

SD= Solid Door
GD= Glass Door



TELR-130-1-SD

TELR-240-1-SD

TELR-250-1-SD
TELR-250-1-GD

TELR-500-1-SD
TELR-500-1-GD

FEATURES

Thermoline manufactures our Economy Style Laboratory Refrigerators to meet the requirements of general purpose laboratory needs, suited for the storage and testing of non critical products where alarms and monitoring are generally not required. We have four models in the economy range from the under bench 140Lt to the upright 500Lt.

- White plastic internal liner featuring CFC free Urethane injection insulation
- White colorbond steel external finish
- Fan forced air supply
- Glass or solid access door (GD or SD)
- Max/min temperature recall feature
- Self closing door (250 litre and 500 litre models)
- Temperature control range between +2°C and +10°C
- 13mm cable monitoring port access
- Quiet operation using hermetically sealed compressor
- Open wire plastic coated shelves
- Adjustable and removable shelves
- Can be positioned up against wall
- Door locks on 130 and 240 litre models

Optional Accessories

Extra shelves
Door locks on 250 and 500 litre models
BMS system for external monitoring
Low temperature safety cut out

Quick Specifications

- Temperature control range +2°C to +10°C
- Fan forced air circulation
- Digital temperature indication
- Maximum/Minimum data recall
- Internal light user selectable on-off

ORDER CODE	DIM. INT HxWxD cm	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg
TELR-130-1-SD	70x46x43	130	84x54x53	2	1A	50
TELR-240-1-SD	127x45x43	240	151x55x56	5	1A	80
TELR-250-1-GD/SD*	105x47x42	220	169x59x62	3	1A	95
TELR-500-1-GD/SD*	142x51x51	500	207x61x77	5	5A	100

(*) When ordering quote **SD** for Solid Door or **GD** for Glass Door.