

**Shakers, Platelet**



TPS-16 (Shaker)

**FEATURES**

The reciprocating shaker is available specifically for platelet bags. It is suitable for the new larger bags now being adopted. The model TPS-16 holds 16 large bags. The unit is attractively finished in a powder coated finish, as is the slide out shelves. It operates at a fixed speed of 60 strokes per minute.

ORDER CODE	TPS-16
CAPACITY	16 large bags or 48 small bags
EXT DIMS (HxWxD cm)	42x52x38
SHELVES	7 slideout and 1 fixed on top
SPEED	60 strokes per minute
ELECTRICAL	240V 50Hz 1A
WEIGHT (kg)	22

**Shaking Incubator, Platelet (+20°C to +24°C)**



TPSRI-16 (Incubator)

**FEATURES**

The bench mounted refrigerated incubator has been designed to suit the TPS-16 Platelet Shaker. It is refrigerated and can operate between 18°C and 24°C in a 30°C ambient. This model is fitted with 150mm dia circular chart recorder (with battery backup) plus Hi Lo alarms and BMS.

The shaker plugs into a power outlet located on the outside of the cabinet via an access port at the rear of the incubator. Forced air circulation provides good air flow in and around the platelet bags.

ORDER CODE	TPSRI-16
CAPACITY	1 TPS-16 shaker
TEMPERATURE RANGE	18°C to 22°C ±1°C
TEMPERATURE CONTROL	Electronic digital indication
ALARMS	Audible deviation alarm Hi and Lo
RECORDER	150mm circular chart (battery backup) with 100 charts
EXT DIMS (HxWxD cm)	85x66x60
INT DIMS (HxWxD cm)	46.5x55x38.5
REFRIGERATION	Air cooled sealed system (auto defrost and condensate disposable)
ELECTRICAL	240V 50Hz 1A
WEIGHT (kg)	60