

## Slide Hot Air Dryer

AUSTRALIAN MADE



THAD-11

### FEATURES

- A small Hot Air Dryer providing a fast efficient method for the drying of small glassware, slides, etc.
- Finished internally in stainless steel with a stainless steel perforated hinged lid.
- The exterior is finished in a PVC coated steel.
- Thermostat control range is +30 to +75°C. Air flow is fan assisted to quickly remove moisture from the enclosure.
- Removable thermometer dial indicator is supplied.

ORDER CODE	DIM. INT LxWxD cm	LITRES CAPACITY	DIM. EXT LxWxD cm	ELECTRICAL 240V 50Hz	PACKED WEIGHT Kg
THAD-11	30x26x14	11	49x30x22	1500W	12