

**Warming Tray (+5°C to +60°C)**

AUSTRALIAN MADE



TWT

**FEATURES**

The Thermoline Warming tray is ideal for the drying of small items such as microscope slides. The large black anodised aluminium warming surface is gently heated by a silicone element attached to the underside. The temperature is set at the digital solid state control.

TEMPERATURE RANGE:      Amb. +5 to +60°C.  
CONTROL ACCURACY:       ±0.25°C

ORDER CODE	HEATED AREA LxW cm	DIM. EXT LxWxH cm	PACKED WEIGHT Kg	ELECTRICAL 240V 50Hz
TWT-D	32x24	43x30x8	3	65W