



GMC-360/361

This compact microcentrifuge is ideal for microfiltration and quick spin downs of various types of microcentrifuge tubes.

The units are supplied with rotors and adaptors to accommodate 2.0ml, 1.5ml, 0.5ml tubes and 0.2ml tube strips.

They are quiet when operating and very easy to use - when turned on, start and stop simply by closing and opening lid.

Model	GMC-361 (with Timer) / GMC-360 (without timer)
Rotor Type	6 place 2.0ml, 1.5ml tube rotor, 6x0.5ml tube adaptor for 2.0ml rotor, 6x0.5ml tube adaptor for 2.0ml rotor
Max Speed	6600rpm for 2.0ml rotor 6300 for 0.2ml rotor
Max RCF	2.000xg for 2.0ml rotor 2200xg for 0.2ml rotor
Lid Material	Polycarbonate
Body/Rotor Material	ABS Resin
Rubber Bottom	Anti-Slip rubber bottom
Dims	160Wx160Dx100Hmm
Power	240V 50Hz