



FR-400-SP

These two spark proof refrigerator models have been specifically designed to allow for the safe storage of volatile samples which are unsafe for storage in standard refrigerators. Both cabinets are completely spark free on the inside only. This ensures that no electrical components, which may cause a spark in the event of a malfunction, are fitted within the workspace.

Both models are constructed internally and externally from stainless steel and are fitted with adjustable open wire shelves for easy storage of samples. Key locks are fitted to the doors to keep any samples safe should the refrigerator be left unattended at any time. Temperature is controlled by the latest controller technology with backlit alphanumeric LCD display. The temperature is also adjustable between 2°C and 15°C. Both models are fitted with heavy duty castors for easy moving.

Optional features include DMLP digital monitor and BMS contacts to connect the refrigerators to a building management system.

NOTE: Only the interior workspace is spark-proof

| Model | FR-400-SP | FR-700-SP |
|---------------------|---|--------------------|
| External Dimensions | 600wx595dx1890h mm | 720wx800dx2050h mm |
| Temperature Range | +2°C to +15°C | |
| Capacity | 400 litres | 700 litres |
| Doors | Single Door - Available in Glass or Solid | |
| Shelves | 5 | 6 |
| Key Locks | Yes | |
| Wheels Kit | Yes | |
| Evaporation Tray | Yes | |
| Weight | 150kg | 240kg |