



MD-01N and MD-02N

Two dry block baths are available to hold single and double dry blocks. Both models include a PI digital controller for easy temperature selection, rapid heat up rates, excellent temperature uniformity of $\pm 0.2^{\circ}\text{C}$ and excellent temperature accuracy of $\pm 0.2^{\circ}\text{C}$. The required temperature can be increased in 0.1°C increments from ambient to 150°C . Both models include a temperature calibration function and a timer from 0 to 9999 minutes. A handy screw in lift out device is supplied for the easy removal of dry blocks.

Block dimensions: 105L x 80W x 50D mm.

Models	MD-01N	MD-02N
No. of Blocks	1	2
Temp. Range	Ambient to 150°C	
Dimensions	260Lx210Wx80D mm	

Block Models	Tube Capacity / Dia	No. of Wells	Depth of Wells
AB-15	15ml / 17mm	12	40mm
AB-10	10ml / 16mm	12	
AB-2.0	2.0ml	20	35mm
AB-1.5	1.5ml	20	34mm parallel
AB-0.5	0.5ml	20	28mm tapered
AB-0.2	0.2ml	30	16mm tapered
AB-MP	96 Well Microtitre Plate (2 blocks)		

Custom blocks available on request.