



ECO-300

The ECO Superfreezer is an economy upright freezer designed to operate between -16 and -24 degrees. This freezer has a digital adjustable controller with audible high temperature and door open alarms, plus an advanced automatic defrosting cycle which prevents ice build up within the chamber, with minimal loss of conditions.

Standard features include four pull out draws and two top fast freeze compartments, reversible front door position and front castors for easy moving.

The Fiocchetti ECO-300 Superfreezer provides constant high performance with minimal power consumption.

Model	ECO-300
Capacity	300 litres
Temperature Range	-16°C ~ -24°C
External Dims (dxwxh mm)	600x630x1600
Reversible Door	Standard
Digital LCD Controller	Standard
Over Temperature Alarm	-12 Degrees
Shelves / Drawers	2 shelve compartments / 4 pull out drawers
Fast Freeze Compartments	2
Energy Consumption	0.73 Kw over 24 hours
Style	Upright
Weight	60 kg
Power	220V 50Hz