



Labor 400 2T
(4002PC)



Labor 400 2T
(DMO0400)
(TGA Approved)

These two temperature refrigerator freezer cabinets are ideal for work areas with limited space available. Both Labor and Vision models have two independently controlled chambers - top section controllable from 0 to +15 degrees and the bottom from -15 to -20 degrees. Both sections have front opening doors which are available in both solid and glass for the Vision model, and top only glass or solid in the Labor model.

Both Labor and Vision models feature the latest controller technology with backlit alphanumeric LCD display and many unique functions including an advanced defrosting feature that allows efficient defrosting cycles which reduces temperature variations and power consumption.

These cabinets are constructed internally and externally from white colorbond steel. All models have internal rounded back corners and a dished bottom. Heavy duty castors are used to keep the cabinets stable when they are being moved.

Vision Model - TGA Approved

The DMD0400 model is approved by the TGA to be used as a medical device in hospitals and medical centres throughout Australia. This model has two 200 litre sections with all the same features which the Labor models offer, however it also includes the DMLP monitor, memory card feature and printer as standard.

Two Temperature Chamber Quick Points

- Slide out aluminium drawers with height adjustment for the bottom section (all models)
- Adjustable slide out shelves for the top section (all models)
- Bottom solid door only and top glass or solid door (Labor models)
- Top and bottom double glazed glass doors or solid doors (Vision model only)
- Lockable self closing door/s (all models)
- Internal lights (all models)
- Automatic condensate disposal (all models)
- Lockable castors (all models)
- DMLP controller with memory card feature and printer (Vision model only)
- BMS contacts to connect to a Building Management System (all models)

Model	2802PC	4002PC	DMO0400
External Dimensions	520wx665x1650 mm	600wx630x1950 mm	600wx655x1950 mm
Temperature Range	+2°C to +15°C - Top -15°C to -20°C - Bottom		
Capacity	140 litres - Top 140 litres - Bottom	200 litres - Top 200 litres - Bottom	
Doors	Upper and Lower Front Opening Doors Labor Model - Bottom Solid Door and Glass/Solid Top Door Vision Model - Top and Bottom Solid/Glass Doors		
Shelves	3 Slide out shelves - Top 3 Slide out drawers - Bottom		
Key Locks	Yes		
Wheels Kit	Yes		
Evaporation Tray	Yes		
DMLP Controller	No	No	Yes
DMLP Printer	No	No	Yes
BMS Connection	Yes	Yes	Yes
Weight	160kg	180kg	180kg