



TS 610

Two large capacity Orbital Shakers available. Both are well engineered and have heavy duty orbital mechanisms that provide smooth operation when under maximum loading and speed.

The top plate is 6mm aluminium and has multiple tapped holes to accommodate flask holders from 50ml up to the heavy flasks sizes: 1 Litre to 6 Litre. Both models have analogue speed display and are attractively finished in a durable grey powder coating.

ORDER CODE	TS 610	TS 620
SHAKING PLATE STYLE	Plate with multiple drilled holes	
SHAKING SYSTEM	Orbital	
RATE OF SHAKING	20-250 RPM	
SHAKING WIDTH	25mm	
SHAKING PLATE (WxD)	600x480	960x600
OUTSIDE DIMS. (WxDxH mm)	600x510x210	960x600x210
	240V 50Hz 2A	
WEIGHT (kg)	52	75

ORDER CODE	TFH-50	TFH-125	TFH-250	TFH-500	TFH-1L	TFH-2L	TFH-3L	TFH-4L	TFH-5L	TFH-6L
FLASK HOLDERS	50 ml	125 ml	250 ml	500 ml	1 litre	2 litre	3 litre	4 litre	5 litre	6 litre
TS-610	80	48	30	20	12	6	6	4	2	2
TS-620	160	96	60	40	24	15	8	8	6	6