



LOB-506D

Three digital oil bath models available, each with a precision digital PID controller, a touch sensitive membrane type front mount controller with timer, alarm and auto-tuning function. Baths are made from stainless steel with a powder coated steel case. These baths are supplied standard with a stainless steel flat lid.

Models	LOB-506D	LOB-511D	LOB-522D
Bath Dims (WxDxH)	300x155x150 mm	300x220x170 mm	400x300x170 mm
Overall Dims (WxDxH)	338x195x400 mm	338x280x400 mm	536x335x400 mm
Capacity	6 litres	11 litres	22 litres
Temp Range	Ambient plus 5°C to 300°C		
Temp Accuracy	±1.5°C		
Controller	Digital PID controller		
Safety	Over temp cut off		
Power	220V 50Hz		



EWK-100

These mini baths feature microprocessor based temperature control with ON/OFF switch and proportional control settings for more precise temperature control. Large dial knob which lets users set the temperature easily. These mini baths can be used in combination with a magnetic stirrer, or as a rotary evaporator or simply as a heater.

Models	EWK-100	EOK-200
Bath Depth	120mm	
Overall Dims (WxDxH)	208 dia x 120 depth	
Capacity	4.5 litres	
Temp Range	ambient +5 to 95°C	ambient +5 to 180°C
Temp Control	Microprocessor control	
Temp Sensor	Type K Thermocouple	
Heater Wattage	750	
Power	220V 50Hz	