



With this colony counter, counting can be achieved by simply touching the petri dish with the counting pen or by pressing the red button located on the front of the counter. The counter confirms each mark with an audible signal and by displaying the counted amount on the LED display. Supplied with counting pen and magnifying glass.

Model	DCC-560
Readout	4 digit LED display
Dish Holder	Petri dish up to 15cm
Illumination	Ring shaped lamp
Count Control	1. manual push button 2. needle point probe 3. marking pen
Power	220V 50Hz
Dimensions	275Wx335Dx175H mm
Weight	5kg