



SHR-2D
(Shown with optional flask tray)

These two reciprocating shakers are ideal for general applications where smooth mixing or intensive mixing is required. Two platform sizes are available, 230x230mm and 350x350mm, both with adjustable speed control by turn dial. Both models have a LCD display which indicates timer duration, current RPM and a power bar graph. Audible alarms to notify the user in the event of an error or timer end.

- Product features include:
- High precision digital speed control.
 - Convenient setting of timer and speed using jog shuttle switch.
 - Storage function for RPM and Timer (99hr 59 min).
 - Locking mode.
 - Smooth start and change of speed.
 - Over load protection.
 - Adjustable shaking width.



Universal Platform



Flask Holders

Universal Platform and Flask Holders

Model	SHR-1D	SHR-2D
Motor	Plate Type BLDC Motor	
Base Platform Size	230x230	350x350
Max Load (kg)	8	10
Distance of Movement (mm)	10, 15, 20 available	20, 30, 40 available
Speed Range (mm)	30 ~ 300	
Speed Control	Programmable operation fully digital feedback control using microprocessor	
Display	Digital LCD Display of RPM timer and Power Bar Graph	
Overload Protection	Automatic stop and smooth speed change	
Dims (wxdxh mm)	250x330x150	400x480x160
Weight (kg)	14	20
Voltage	220V 50Hz	

Model	Description
SHO-2D-SP100	Universal platform for SHR-1D, 230wx230d mm
SHO-2D-SP110	Rubber Mat platform for SHR-1D, 270wx270d mm
SHO-2D-SP200	Universal platform for SHR-2D, 250wx250d mm
SHO-2D-SP210	Rubber mat platform for SHR-2D, 370wx370d mm
TFH-50	Flask holder for shaker, 50ml
TFH-100	Flask holder for shaker, 100ml
TFH-250	Flask holder for shaker, 250ml
TFH-500	Flask holder for shaker, 500ml