



WOV-30

Two digital control vacuum ovens are available in 30 and 70 litre capacities. Both models are suitable for drying, baking, conditioning, vacuum embedding, moisture testing and age testing. All controls are conveniently located on the front panel for easy access. These models have a robust design and provide reliable use over many years of service.

- Electronic digital control technology for accurate control
- User select auto or manual run applications
- Large Vacuum gauge for easy viewing
- High temperature rating silicon door gasket
- Safety glass view window
- Storage ability for alarms, timer on digital control panel
- Lock mode
- Stainless steel 304 internal sealed shell
- Removable inner shelves

Safety features include over temperature protection, over current protection and sensor fail indication.

Model	WOV-30	WOV-70
Capacity	30 litres	70 litres
Internal Size (wxdxh mm)	300x330x300	400x435x400
External Size (wxdxh mm)	470x455x575	570x560x675
Heater Capacity	1000w x 2	550w x 4
Airflow	Natural Convection	
Temp. Range	Ambient to 200°C	
Temp. Accuracy	±2% of the current temperature	
Temp. Controller	Digital Fuzzy Controller	
Vacuum Range	0 to 750 mmHg	
Timer Range	Digital 99 Hours to 59 minutes	
Shelves	Two removable compartments	Three removable compartments
Material Internal	Stainless Steel	
Material External	Powder coated steel	
Insulation	Ceramic Wool	
Door Packing	High Temperature Silicone Rubber	
Window	High Temperature Safety Glass	
Safety Device	Over temperature, over current protection	
Power	AC 230V, 50Hz	