



RK-1D

These two rockers are ideal for general applications where smooth mixing is required. Two platform sizes are available, 290x200mm and 300x300mm, both with adjustable speed control by turn dial. Both models have a LCD display which indicates timer duration, current RPM and a power bar graph. Audible alarms to notify the user in the event of an error or timer end.

Product features include:

High precision digital speed control.

Programmable operation modes supported 10 steps , 6 programs.

Convenient setting of timer and speed using jog shuttle switch.

Storage function for RPM and Timer (99hr 59 min).

Locking mode.

Smooth start and change of speed.

Over load protection.

Model	RK-1D	RK-2D
Platform Dims (wxd mm)	290x200	300x300
Overall Dims (wxdxh mm)	290x240x150	300x330x150
Tilt Angle	9° (default)	
Motor	35W geared DC motor	
Speed Range	5 - 50 digital 100 Hr digital	
Timer	100 Hr digital	
Display	Digital LCD display	
Safety	Over load detector	
Options	Optional platform with 4 pillar	
Power	220V 50Hz	
Alarms	Emergency alarm (buzzer)	