

**Baths, Boiling (+100°C)**

AUSTRALIAN MADE



TBB-8

**FEATURES**

This boiling bath is designed to operate at +100°C and is fitted with removable concentric rings ranging in size from 30, 50, 75 and 100mm in diameter to support round based flasks of varying sizes.

- Seamless stainless steel bath with fibreglass insulation and high capacity calrod element for rapid heat up.
- Temperature control is by an adjustable Simmerstat Energy Regulator.
- Fitted with a constant water level device and can be connected to an adjacent tap.
- The bath needs to be able to drain constantly to maintain the correct water level.
- This model has a Manual Reset Safety Thermostat.

**Safety Devices**

- ✓ Constant water level device
- ✓ Manual Reset Safety Thermostat fitted as standard

ORDER CODE	DIM. EXT LXWXD	LITRES CAPACITY	Nº OF RINGS	ELECTRICAL 240V 50HZ	APPROX WEIGHT Kg
TBB-8	62x33x13	9	8	2000W	10