

**Baths, Circulated (+100°C with Heater Circulator)**

AUSTRALIAN MADE



TWBC-24 shown with optional TU3

Note: Heater Circulator pictured is optional with the TWBC-Series

ORDER CODE	DIM. INT USEABLE LxWxD	LITRES CAPACITY	DIM. EXT LxWxD	APPROX WEIGHT	OPTIONAL GABLE LIDS	OPTIONAL FLAT LIDS
TWBC-12	16x23x20cm	12	33x27x23cm	6kg	PCLC-12 (polycarbonate)	SSFLC-12 (stainless steel)
TWBC-24	36x29x20cm	24	53x33x23cm	8kg	PCLC-24 (polycarbonate)	SSFLC-24 (stainless steel)
TWBC-48	48x49x20cm	48	65x53x23cm	15kg	SSGLC-48 (stainless steel)	SSFLC-48 (stainless steel)

**FEATURES**

Thermoline have a quality range of general purpose circulated water baths. These baths are Australian designed for safety, reliability and longevity.

- Insulated Circulated Bath standard sizes 12, 24, and 48 litre capacity
- Seamless Stainless Steel tank
- Radial corners (quick and easy to clean or disinfect)
- Carry handles fitted as standard
- Custom baths can be ordered to suit your application
- Temperature range when using TU-series of Heater Circulators (see page B-4) ambient +5°C to +100°C
- Extended temperature range of -30°C to +80°C circulated baths available



TWBC-12 shown with optional TU3

- Optional Accessories**
- ✓ Polycarbonate gable lids to suit TWBC-12 and TWBC-24
  - ✓ Stainless steel gable lid to suit TWBC-48
  - ✓ Stainless steel flat lids available to suit all models
  - ✓ TU-Series of heater circulators to suit all your application