

Shaking Thaw Bath (Ambient +5°C to +45°C)

AUSTRALIAN MADE



Raised rim to reduce splashing



Special thaw tray

FEATURES

- Specially designed by Thermoline for thawing frozen products
- Large capacity to accommodate 12 x500ml frozen plasma bags
- Splash protection rim
- Fast Thawing time reduced to 12 minutes on the adjustable timer (full load 30 minutes)
- Digital precision of control temperature better than $\pm 0.2^{\circ}\text{C}$
- Dual digital display for quick and easy viewing of set temperature and bath actual temperature
- Timer selection user adjustable from 1 min to 999 minutes with audible time out function
- Unique thaw tray allows for super quiet operation and perforated sides, allow the middle bags to be thawed at the same time
- Locking springs keeps all bags in the correct position

Safety Devices

- ✓ Over temperature cut out set at +45°C
- ✓ Auto over temperature safety (Failsafe)
- ✓ Audible time out signal
- ✓ Control error codes
- ✓ Manual Reset Safety Device (MRST)

TEMPERATURE RANGE:
Ambient +5°C to +45°C
(Recommended operating temperature 37°C)

CONTROL STABILITY: $\pm 0.25^{\circ}\text{C}$

ORDER CODE	DIMS EXT (HxWxD CM)	TRAY (LxWxH cm)	BAG CAPACITY	ELECTRICAL 240V 50Hz	WEIGHT kg
TSBT-20	66x33x29	37x27x5	12x500ml	2100W	22