

Freezers, Economy (-15°C to -20°C)

AUSTRALIAN MADE



TEUF-181



TEUF-90

FEATURES

These two economy style upright freezers are ideal for samples requiring storage at -20°C. Both models are fitted with a digital LAE controller, which has high and low alarms, and are both supplied standard with manual defrost (manual defrost ensures the temperature remains stable during operation). Shelves are in a fixed position and clear plastic storage containers are supplied.

- Key lock fitted to door.
- Added benefit of digital control with indication of internal conditions.
- Easy installation - simply plug in and set temperature.

Optional Features

POWER GUARD, offering protection against low voltage, high voltage spikes and limited lighting

| Model | TEUF-90 | TEUF-181 |
|-----------------------|---|-----------------|
| Capacity | 90 Litres | 181 litres |
| Temp. Range | -15°C to -20°C | |
| Temp. Uniformity | ±0.2 | |
| Storage Containers | 4 clear plastic | 6 clear plastic |
| Ext. Dims. (hwxwd mm) | 870x540x600 | 1440x545x570 |
| Int. Dims. (hwxwd mm) | 650 / 430x400x400 | 1240x405x380 |
| External Material | White enamel | |
| Internal Material | White PS | |
| Electrical | 240v 50Hz, Max 4A | |
| Weight | 40kg | 60 kg |
| Options | External connection for BMS (in-house monitoring) | |