

Upright Freezers, Superartic (-20°C to -40°C with auto defrost)



Superartic 250



Superartic 600

FEATURES

These two new upright freezers in the Superartic series have been specially designed for products which need to be stored at ultra low temperatures. Both freezers have an adjustable temperature range from -20°C to -40°C.

The cabinets are constructed internally with stainless steel and externally from white colorbond steel. Polyurethane foam insulation is used to help maintain the internal temperature. Heavy duty castors are used to keep the freezers stable when they are being moved. Cabinets constructed with stainless steel exterior are available at an extra cost.

Both models feature the latest controller technology with backlit alphanumeric LCD display for easy temperature adjustments and a password protect on/off option which prevents set point modifications without authorisation. Audible and visual alarms including recording of high/low temperature, power failure, door open and dirty condensor alarms notify the user should a problem occur during operation.

Product Quick Points

- Reversible front door to suit your location
- Lockable doors
- Adjustable shelves or optional drawers
- Forced air circulation for quick freezing
- Advanced automatic defrost and condensate disposal
- BMS contacts to connect to a Building Management System

Safety Features

- ✓ Alarms fitted standard
- ✓ Over current protection
- ✓ Lockable castors

Optional Features

Circular Chart Recorder

ORDER CODE	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 230V 50Hz	APPROX WEIGHT Kg
Superartic 250	250	158x67x83	4	5A	150
Superartic 600	600	208x72x76	7	5A	280