

Economy Incubators (Ambient +5°C to 60°C)

AUSTRALIAN MADE



FEATURES

Thermoline have four economy incubators available, all designed and manufactured in house. These models are economical and very easy to use making them an ideal choice for all non critical samples requiring incubation. A wide range of fan forced models are also available, should tighter control conditions and larger volumes be required (see page 14).

Standard Features includes

- Stainless steel internal liner (easy to clean and non corrosive)
- Removable open wire chrome plated shelves
- Magnetic door gasket means positive door closure
- Light weight and portable
- Digital Microprocessor control complete with dual digital display
- Control error indication
- Unique “Thermowall” heating system with electronic switching
- Control range ambient +5°C to +60°C
- Control Accuracy at set point ±0.25°C or better

Safety Devices

- ✓ Auto overtemperature safety (Failsafe)
- ✓ Electrical over current protection
- ✓ Control error codes

Optional Accessories

- IGD- Internal glass door installed (prevents loss of internal conditions when access door is opened)*
- GVW- Glass View Window installed in access door*

ORDER CODE	DIM. INT HxWxD cm	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg	OPTIONAL INNER GLASS DOOR	OPTIONAL GLASS VIEW WINDOW
TEI-13G	22x25x24	13	37x30x30	2	45W	8	IGD	GVW
TEI-43G	50x29x29	43	65x34x35	4	90W	16	IGD	GVW
TEI-70G	60x29x40	70	75x34x46	5	130W	24	IGD	GVW
TEI-100G	65x39x40	100	83x46x47	6	180W	32	IGD	GVW