

## Economy Refrigerated Incubators (+5°C to +40°C / +60°C)

AUSTRALIAN MADE



TREI-150-1-VW



TREI-220-1-SD



TREI-490-1-VW

### FEATURES

Thermoline has three economy refrigerated incubators available in 150, 220 and 490 litres in capacity. All models feature precision temperature control with fan forced air circulation. Adjustable temperature range available for the 150 litre and 490 litre models is from +5°C to +60°C, with the cabinet capable of running at 0°C for a maximum of three days before defrosting is required. The 220 litre model has an adjustable temperature range limited to +5°C to +40°C.

Both 150 litre and 490 litre models have a powder coated steel exterior with a polished stainless steel interior. The 220 litre model has a durable white plastic interior and white enamelled steel exterior.

- Self contained air cooled refrigeration system offering quiet operation and ozone friendly refrigerant
- Visual indication of set point & internal temperatures
- Precise temperature control  $\pm 0.20^{\circ}\text{C}$  at set point with  $0.1^{\circ}\text{C}$  readability
- Cable port hole
- Adjustable open wire chrome plated shelves
- Over temperature safety cut out fully adjustable within range
- Sliding window cover for 150 and 490 litre models.

#### Safety Devices:

- ✓ CFC free refrigeration
- ✓ Auto overtemperature safety (Failsafe)
- ✓ Over current protection

#### Available temperature ranges

- |                      |  |
|----------------------|--|
| <b>TREI-150-1-VW</b> | +5°C to +60°C (0°C for a limited period) |
| <b>TREI-220-1-SD</b> | +5°C to +40°C                            |
| <b>TREI-490-1-VW</b> | +5°C to +60°C (0°C for a limited period) |

#### Available Options:

*POWER GUARD, offering protection against Low Voltage, High Voltage Spikes*  
*Alarms with Audible & Visual indication & remote contacts*  
*Auto dialler on alarm action*  
*Extra shelves*  
*Additional cable entry points*

ORDER CODE	DIM. INT HxWxD cm	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg
TREI-150-1-VW	50x60x50	150	89x68x71	3	5A	90
TREI-220-1-SD	127x44x40	220	155x54x55	4	4A	95
TREI-490-1-VW	120x70x58	490	177x79x80	6	8A	188