

Laboratory Refrigerators, Economy (+4°C to +10°C)

AUSTRALIAN MADE

Standard features include:

- Internal Forced Air Flow
- Digital Temperature Display with Alarms
- Max/Min Temperature Logging
- Internal Light
- Cable Port
- Door Key lock
- Frost Free Operation
- CFC Free Refrigerant
- Adjustable/Removable Shelves



TEL R-130-1-SD



TEL R-440-1-GD



TEL R-740-2-GD

FEATURES

Five Economy Laboratory Refrigerators are available up to 740 litres in capacity to meet the general purpose laboratory requirements for storing and testing of non critical products. A further four models are available up to 1150 litres capacity for increased storage requirements. (see pages R-2)

Temperature is precisely controlled electronically between +4°C and +10°C in ambients up to 32°C. The digital display incorporates high and low temperature alarms, which emit visual and audible signals to warn the operator of an alarm condition. Connection to a Building Management System (BMS) is available as an optional extra.

A Low Temperature Protective Cut Out (LTPC) can be fitted as an optional extra to prevent products inside the refrigerator freezing. This independent device will turn off the refrigerator system should the control fail in the cooling mode.

Models TEL R-130, TEL R-240 and TEL R-270 have a white metal exterior and a white moulded plastic interior with curved corners. Larger models TEL R-440 and TEL R-740 have a grey metal exterior and interior complete with stainless steel floor and white Colorbond walls. Insulation is high efficiency injected polyurethane foam.

All models are supplied standard with self closing key lockable doors. Models TEL R-130, TEL R-240 and TEL R-270 have insulated solid doors whilst the larger models TEL R-440 and TEL R-740 have double glazed glass doors which enable the products to be viewed without any loss of conditions inside the refrigerator.

OPTIONS

- LTPC (Low temperature protective cut out)*
- Additional Shelving*
- BMS (Voltage free contacts allow connection for remote monitoring of power loss, low and high temp. alarms)*
- Phone Dialler (TPD-1) Analogue line only.*
- POWER GUARD offering protection against Low Voltage, High Voltage Spikes*

ORDER CODE	DIM. INT HxWxD cm	LITRES CAPACITY	DIM. EXT HxWxD cm	SHELVES	ELECTRICAL 240V 50Hz	APPROX WEIGHT Kg
TEL R-130-1-SD	70x46x43	130	84x54x53	2	1A	50
TEL R-240-1-SD	127x45x43	240	151x55x56	5	1A	80
TEL R-270-1-SD	146x46x43	270	165x53x57	5 / 1 drawer	1A	90
TEL R-440-1-GD	130x50x52	440	200x65x75	4	4.2A	150
TEL R-740-2-GD	130x85x52	740	200x100x75	4	4.2A	240