

Warming Tray (+5°C to +60°C)

AUSTRALIAN MADE



TWT

FEATURES

The Thermoline Warming tray is ideal for the drying of small items such as microscope slides. The large black anodised aluminium warming surface is gently heated by a silicone element attached to the underside. The temperature is set at the digital solid state control.

TEMPERATURE RANGE:	Amb. +5 to +60°C.
CONTROL ACCURACY:	±0.25°C

ORDER CODE	HEATED AREA LxW cm	DIM. EXT LxWxH cm	PACKED WEIGHT Kg	ELECTRICAL 240V 50Hz
TWT-D	32x24	43x30x8	3	65W