

AUTOMATIC TUBE STRIPPER

(HAND HELD)

Image: TS08



Automatic Tube Stripper

- Designed for left hand or right hand use
- Fast and complete stripping of blood tubes
- Ideal for all commercial blood bag systems on the market
- Tubing automatically centered through rollers with guide
- Changeable stripping speed
- High efficiency and low energy consumption



Rollers with tubing

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au

 **Thermoline**
S C I E N T I F I C
H I G H E R S T A N D A R D S I N S C I E N C E

PRODUCT DETAILS

This automatic blood tube stripper allows users to quickly and easily remove left over blood within blood tubes by squizing the tube between two rollers. When the lever is pressed, the rollers turn forcing the excess blood through the tube and out the end.

SPECIFICATIONS

Model	TS08
Overall Dimensions (handle open mm)	51x79x175
Maximum Tube Diameter (mm)	145x79x175
Stripping Wheel Diameter (mm)	6.5
Stripping Wheels Angle Opening	25°
Tube Stripper Body Diameter (mm)	40
Maximum Stripping Speed (mm / sec)	400
Weight (g)	600
DC Motor Power	5-24 V / 1 V Step

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au