

# CIRCULATED WATER BATHS

(AMBIENT +5°C TO +100°C)

Image: TWBC-24 shown with optional TU3



## Thermoline Digital Circulated Water Baths

- 12, 24 and 48 litre models available
- Bath constructed of high grade stainless steel
- Insulated walls to reduce heat loss
- Seamless stainless steel tank with radial corners
- Carry handles for easy moving
- Custom made baths available on request
- Australian made

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email [info@thermoline.com.au](mailto:info@thermoline.com.au)

Web [www.thermoline.com.au](http://www.thermoline.com.au)

 **Thermoline**  
S C I E N T I F I C  
HIGHER STANDARDS IN SCIENCE

# PRODUCT DETAILS

Thermoline's quality range of stainless steel laboratory water baths available in 12, 24 and 48 litres capacity. All baths are constructed of stainless steel, and are insulated to reduce heat loss. Accurate temperature control is available when used together with the TU series heater circulators.

Optional Accessories
Polycarbonate gable lids to suit 12 and 24 litre models
Stainless steel gable lid to suit 48 litre model
Stainless steel flat lids available to suit all models
TU Series heater circulators to allow heating from ambient +5°C to 100°C

- For more information on Thermoline's range of Unistat Heater Circulators, please refer to our "Digital Heater Circulators" page.
- Information on Thermoline's Uncirculated Baths can be found on our "Digital Uncirculated Baths" page.

## SPECIFICATIONS

Model	TWBC-12	TWBC-24	TWBC-48
Capacity	12 Litres	24 Litres	48 Litres
Internal Dimensions (WxDxH mm)	160x230x200	370x300x200	490x500x200
External Dimensions (WxDxH mm)	330x270x230	530x330x230	650x530x230
Optional Gable Lids	PCLC-12 (polycarb)	PCLC-24 (polycarb)	SSGLC-48 (stainless)
Optional Flat Lids	SSFLC-12 (stainless)	SSFLC-24 (stainless)	SSFLC-48 (stainless)
Weight (kg)	6	8	15

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

**Head Office Phone** (02) 9604 3911  
**Head Office Fax** (02) 9725 1706

**Email** info@thermoline.com.au  
**Web** www.thermoline.com.au

# CIRCULATED WATER BATHS

(AMBIENT +5°C TO +100°C)

Image: WCB-11 with flat lid



## **Daihan Digital Circulated Water Baths**

- 6, 11 and 22 litre models available
- Bath constructed of high grade stainless steel
- Economic and easy to operate
- Seamless stainless steel tank with radial corners
- Flat and gabel lids available for all models

**Head Office Phone** (02) 9604 3911

**Head Office Fax** (02) 9725 1706

**Email** [info@thermoline.com.au](mailto:info@thermoline.com.au)

**Web** [www.thermoline.com.au](http://www.thermoline.com.au)

# PRODUCT DETAILS

These three WiseCircu circulated water baths are powerful and economic. All models are fitted with a digital fuzzy control system with backlit LCD display and a turn dial for easy temperature and timer adjustments. Pump out facility to allow connection to external instruments/baths is supplied standard with all three models. Temperature compensation is supported for user calibration.

All baths are supplied standard with a flat stainless steel lid.

## Optional Accessories

Gable lids available for all models

# SPECIFICATIONS

Model	WCB-6	WCB-11	WCB-22
Capacity	6	11	22
Bath Dimensions (WxDxH mm)	300x155x150	302x240x150	500x295x150
Overall Dimensions (WxDxH mm)	338x195x400	338x280x400	536x335x400
Temperature Range	Ambient +5 °C to 100 °C		
Controller	Digital Fuzzy Control System		
Safety	Adjustable Over Temperature, Sensor Error and Low Water Level		
Heater (W)	800	1200	2000
Pump	9 Litres / Min		
Timer	99 Hr 59 min		

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

**Head Office Phone** (02) 9604 3911

**Head Office Fax** (02) 9725 1706

**Email** [info@thermoline.com.au](mailto:info@thermoline.com.au)

**Web** [www.thermoline.com.au](http://www.thermoline.com.au)