

COLONY COUNTER

(WITH MAGNIFY GLASS AND MARKING PEN)

Image: DCC-560



Colony Counter with Magnify Glass and Pen

- 4 digital LED display
- Suitable for petri dishes up to 150mm
- Counting can be done via marking pen or manual push button
- Magnify glass and marking pen supplied standard

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au

 **Thermoline**
S C I E N T I F I C
H I G H E R S T A N D A R D S I N S C I E N C E

PRODUCT DETAILS

With this colony counter, counting can be achieved by simply touching the petri dish with the counting pen or by pressing the red button located on the front of the counter. The counter confirms each mark with an audible signal and by displaying the counted amount on the LED display. Supplied with counting pen and magnifying glass.

SPECIFICATIONS

Model	DCC-560
Readout	4 Digital LED display
Dish Holder Size	Petri dish holder up to 150mm
Illumination	Ring shaped lamp
Counting Control	manual push button or marking pen
Power	220V 50Hz
Dimensions (WxDxH mm)	275x335x175
Weight (kg)	5

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au