


DIGITAL ROCKERS

(BENCH TOP)

Image: RK-2D



Daihan Digital Rocker

- Programmable operation modes supported 10 steps , 6 programs
- Convenient setting of timer and speed using jog shuttle switch
- Storage function for RPM and Timer (99hr 59 min)
- Locking mode
- Smooth start and change of speed

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email info@thermoline.com.au

Web www.thermoline.com.au

PRODUCT DETAILS

These two rockers are ideal for general applications where smooth mixing is required. Two platform sizes are available, 290x200mm and 300x300mm, both with adjustable speed control by turn dial. Both models have a LCD display which indicates timer duration, current RPM and a power bar graph. Audible alarms to notify the user in the event of an error or timer end.

Safety Devices

Over current protection

SPECIFICATIONS

Model	RK-1D	RK-2D
Platform Dimensions (WxD mm)	290x200	300x300
Overall Dimensions (WxDxH mm)	290x240x150	300x330x150
Tilt Angle	9°	
Motor	35W geared DC motor	
Speed Range (rpm)	5 to 50	
Timer (hr)	100	
Display	LCD display with backlight	
Safety	Over load protector	
Options	Optional second platform with four support pillars	
Power	220V 50Hz	
Alarms	Buzzer on error	

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email info@thermoline.com.au

Web www.thermoline.com.au