

ECONOMY FREEZER WITH AUTO DEFROST

(-16°C TO -24°C)

Image: ECO-300



Economy Freezer with Automatic Advanced Defrosting

- Advanced defrosting cycle with minimal loss of conditions
- Digital LCD display supplied standard
- Supplied with two shelf compartments and four pull out drawers
- Two fast freezer compartments supplied standard
- Adjustable temperature down to -24°C
- Castors fitted for easy moving

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au



Thermoline
S C I E N T I F I C

HIGHER STANDARDS IN SCIENCE

PRODUCT DETAILS

The ECO Superfreezer is an economy upright freezer designed to operate between -16 and -24 degrees. This freezer has a digital adjustable controller with audible high temperature and door open alarms, plus an advanced automatic defrosting cycle which prevents ice build up within the chamber, with minimal loss of conditions.

Standard features include four pull out draws and two top fast freeze compartments, reversible front door position and front castors for easy moving.

The Fiocchetti ECO-300 Superfreezer provides constant high performance with minimal power consumption.

SPECIFICATIONS

Model	ECO-300
Capacity (litres)	300
Temperature Range	-16°C to -24°C
External Dimensions (DxWxH mm)	600x630x1600
Reversible Door	Yes
Digital LCD Controller	Yes
Over Temperature Alarm	High alarm at -12°C
Shelves / Drawers	2 shelve compartments and 4 pull out drawers
Fast Freezer Compartment	2 compartments
Energy Consumption	0.73 Kw over 24 hours
Style	Upright
Weight (kg)	60
Power	220V 50Hz

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email info@thermoline.com.au

Web www.thermoline.com.au