

LARGE CAPACITY WATER BATHS

(WITHOUT CIRCULATORS)

Image: TLWB-30 with gable lid



Thermoline Large Capacity Water Baths

- Large laboratory and industrial water baths
- Custom baths can be ordered to suit most applications
- Stainless steel inner tank
- Height adjustable perforated floor
- Fully insulated
- Baths to be fitted with optional heater circulator



Optional heater circulator sold separately

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au

 **Thermoline**
S C I E N T I F I C
H I G H E R S T A N D A R D S I N S C I E N C E

PRODUCT DETAILS

Bath Model	Recommended Heater Circulators
TLWB-30	TU-1, TU-3 and TU-4
TLWB-80	TU-5

SPECIFICATIONS

Model	TLWB-30	TLWB-80
Capacity	30 Litres	80 Litres
Internal Dimensions (WxDxH mm)	500x300x200 (useable)	600x350x350 (useable)
External Dimensions (WxDxH mm)	730x360x310	880x460x510
Weight (kg)	15	40

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email info@thermoline.com.au

Web www.thermoline.com.au