

# PREMIUM REFRIGERATED INCUBATORS

(+5°C TO +45°C)

Image: TRI-1175-2-SD



## **Thermoline Premium Refrigerated Incubators**

- Premium custom air handling duct system providing uniform temperatures
- Microprocessor controller displaying internal & set point values of 0.1°C readability
- Premium stability  $\pm 0.2^{\circ}\text{C}$  temporal  $\pm 1.75^{\circ}\text{C}$  or better spatial variation
- Fan forced air circulation
- Premium toughened glass internal viewing doors
- Stainless steel internal liner
- Adjustable open wire chrome plated shelving
- Key door locks supplied
- Heated door fascias which prevent sticking doors

**Head Office Phone** (02) 9604 3911

**Head Office Fax** (02) 9725 1706

**Email** [info@thermoline.com.au](mailto:info@thermoline.com.au)

**Web** [www.thermoline.com.au](http://www.thermoline.com.au)

 **Thermoline**  
S C I E N T I F I C  
H I G H E R S T A N D A R D S I N S C I E N C E

# PRODUCT DETAILS

Thermoline's range of premium refrigerated incubators are ideal for any laboratory that require large capacity cabinets which provide stable and uniform incubation. All premium models have a controllable temperature range between +5°C and 45°C and all have fan forced air circulation via a customised rear air duct system that allows conditioned air to flow across each shelf position. A 13mm port hole is provided to allow sensors or probes to be inserted into the cabinet for monitoring.

## Safety Devices

Over current protection

Controller error codes

Auto over temperature safety (Failsafe)

Key lockable doors

Lockable castors

## Optional Accessories

Remote alarm contacts

POWER GUARD - Offers protection against low voltage, high voltage spikes and limited lightning.

USB data loggers

Additional shelves

# SPECIFICATIONS

Model	TRI-750-2-SD	TRI-1175-2-SD
Internal Dimensions (HxWxD mm)	1290x990x500	1350x1430x550
External Dimensions (HxWxD mm)	1990x1110x730	2010x1550x780
Capacity (litres)	750	1175
Shelves	4	4
Doors	2 Solid Doors	
Electrical (Amps)	8	10
Weight (kg)	180	230

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

**Head Office Phone** (02) 9604 3911

**Head Office Fax** (02) 9725 1706

**Email** [info@thermoline.com.au](mailto:info@thermoline.com.au)

**Web** [www.thermoline.com.au](http://www.thermoline.com.au)