

ROLLER MIXERS

(BENCH TOP)

Image: 210RM



Bench Top Roller Mixers

- Simultaneous rotating and rocking action
- Rollers made from non-slip plastic
- Quiet motor with dural gears for smooth mixing
- Base plate made from stainless steel

Note: New models available with adjustable turn dial speed control from 12 to 70 RPM (9 steps). Quote "C" at the end of the item code when ordering.

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au

PRODUCT DETAILS

Simultaneous rotating and rocking action provides fast and gentle mixing. All rollers are made of special plastic for non-slip and for easy cleaning - 32mm dia x 320mm long rollers. The quiet motor, durable plastic gears, smooth mixing action result in a low noise level. Base plate is made of stainless steel for easy maintenance.

SPECIFICATIONS

Model	205RM*	210RM*
Rollers	5	10
Speed (rpm)	40	50
Roller Length (mm)	325	
Dimensions (LxWxH mm)	480x220x110	480x420x110
Weight (kg)	4.5	8.5

*Quote "C" at the end of item code for adjustable speed model.

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email info@thermoline.com.au

Web www.thermoline.com.au