


TEST TUBE ROTATOR

(DIGITAL)

Image: RT-10



Bench Top Roller Mixers

- Feedback control ensures speed is maintained even when under maximum load
- Smooth adjustable speed range from 0 to 60rpm
- Adjustable angle from 0 to 90 degrees
- Single or double disc option
- Can spin 50ml, 15ml and 1.5ml tubes
- LCD display with 99 hour timer and alarms

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au

PRODUCT DETAILS

The WiseMix RT is a new digital feedback test tube rotator which is small and easy to operate. This rotator has a smooth adjustable speed range from 0 to 60 rpm and an adjustable angle from 0 to 90 degrees. Standard disc sizes are available to spin 50ml, 15ml and 1.5ml conical tubes. Using the screw in bars, this rotator can spin one or two discs together (optional). Features include digital LCD display, 99 hour timer with alarm and locking mode.

Available Discs
WRT100 : 8x50ml tubes
WRT110 : 16x15ml tubes
WRT120 : 32x1.5ml tubes
WRT140 : 4x250ml tubes

SPECIFICATIONS

Model	RT-10
Angle Range	0 to 80 degrees
Speed Range (rpm)	0 to 60
Controller	Digital Feedback with Jog Switch
Display	LCD Display
Timer	99 hour 59 min
Dimensions (WxDxH mm)	200x250x260
Weight (kg)	8.5

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email info@thermoline.com.au

Web www.thermoline.com.au