

VORTEX MIXERS

(ANALOGUE)

Image: 250VM and 260VM



Bench Top Vortex Mixers

- Variable speed control via turn dial
- Suitable for slow speed up to high speed vortexing
- Vortexing speed from 0 to 3000 rpm
- Touch mode and continuous modes of operation
- Single test tube cup head and flat disc style head available
- Light weight and easy to operate

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au

 **Thermoline**
S C I E N T I F I C
H I G H E R S T A N D A R D S I N S C I E N C E