



TRC-500



TRC-1500



TRC-3000

FEATURES

Three Recirculating Cooler models available.

TRC-500 is a unit with a cooling capacity of 500 watts with delivery water temperature of +20°C. It has its own internal pump which is a magnetic drive centrifugal type with a capacity of 20L/min at 1m head. It is suitable to cool water baths or external heat exchanges down to -5°C with the use of a suitable antifreeze. Thermostat control to $\pm 1^\circ\text{C}$ or better.

TRC-1500 is a unit with a cooling capacity of 1500 watts with the delivery water temperature of +20°C. This unit has a powerful bronze centrifugal pump capable of pumping 10L/min up to 8m head. This model is suitable for large baths and heat exchanges with restrictive circuits. Ideal for glass condensers and distillers.

Thermostat control to $\pm 1^\circ\text{C}$ or better.

TRC-3000 is a unit with a cooling capacity of 3000 watts and with a delivery water temperature of +20°C. The pump is a stainless steel centrifugal pump capable of pumping 20L/min up to 20m head. Minimum water temp with this pump is 10°C.

Digital solid state temperature control with accuracy of $\pm 0.5^\circ\text{C}$ or better.

Larger models available on request.

ORDER CODE	REFRIG. UNIT HP.	MAX COOLING CAPACITY WATTS @ (AMB 24°C)			PUMP CAPACITY	EXT DIMS. HxLxD cm	ELECTRICAL	APPROX WEIGHT Kg
		+ 5°C	+ 10°C	+ 20°C				
TRC-500	1/4	380	430	500	20 L/min @ 1M HEAD	27x57x37	5A 240V 50Hz	33
TRC-1500	3/8	500	800	1500	10 L/min @ 8M HEAD	61x57x37	8A 240V 50Hz	55
TRC-3000	7/8	2300	2600	3000	20L/min @ 20M HEAD	100x65x60	10A 240V 50Hz	70

*Specifications are subject to change without notification.

QUEENSLAND SALES

Phone: (07) 3260 5243 | Fax: (07) 3266 3978.
 Email: qldsales@thermoline.com.au

NSW SALES AND MANUFACTURING

Phone: (02) 9604 3911 | Fax: (02) 9725 1706.
 Email: sales@thermoline.com.au

VICTORIA SALES

Phone: (03) 9354 5788 | Fax: (03) 9354 5799
 Email: vicsales@thermoline.com.au