



Cell washer (manual) 12 place Centrifuge 2030. Brushless drive



- LED display of Speed & Time
 - Digital control of parameters
 - Automatic soft start
 - Ten Acceleration & Deceleration rates
 - Ten Memory / Preset Program
 - 1-90 second Pulse (Quick spin)
 - CE Marked (see safety features)
 - Sensors Tacho, Imbalance & Lid
 - Stainless steel bowl & barrier ring
- **Complete with Rotor & Accessories**



Speed 1000-3500 Rpm **Max Rcf** 1400G **Timer** 0-30mins **Size** 23x26x39Cm **Weight** 13.75Kg

2030 Rotors

Adaptors supplied with rotor

12 of each

Tube size 75 x 7mm

Tube size 100 x 12mm



5450

Complete with adaptors fitted



5450

Complete with 24 adaptors

