

# C2006 Fixed angle up to 6000RPM

*"functional simplicity"*

**Speed** 500-6000Rpm

**Max Rcf** 3800G

**Timer** 0-99Minutes

**Dims HWD** 27x33x46Cms

**Weight** 16.5Kilos



## Rotors Available



**BR5412**

max speed 6000 Rpm 3,800G  
**12 x 15ml tube, fixed angle rotor**  
Tube size maximum H120 x W17 mm  
Cap size dependant



**8009**

ADAPTORS for BR5412 (above) 5-10ml  
Tube size maximum H100 x W13 mm



**BR5406**

max speed 6000 Rpm 3,800G  
**6 x 50ml Tube, fixed angle rotor**  
Tube size maximum H115 x W30mm



**6006**

ADAPTORS for BR5406 (above) 25ml  
Tube size maximum H115 x W24 mm

**6007**

ADAPTORS for BR5406 (above) 15ml  
Tube size maximum H120 x W17 mm